



## eXPRESS-LINE Split Tube Furnace

- Quick Delivery
- ID Ranges From 3" To 6.0"
- 12" Long To 36" Long Heated Lengths
- Single Or Three Zone Configuration
- Available With Standard 'QD' Control System
- Includes Touchscreen SmartControl



## Solid Tube Furnace

- 3" ID X 36" Long Heated Chamber
- Three Zone Operation
- Temperatures up to 1250 °C (2282 °F)

## Hinged Split Tube Furnace

- Silicon Carbide Elements
- Temperatures up to 1500 °C (2732 °F)
- Vertical or Horizontal Operation



## Vertical 3 Zone Split Tube Furnace With Stand

- Hinged Furnace Halves
- Heated Chamber: 8" ID X 72" High
- Optional Pneumatic Opening/Closing Mechanism



## Rotary Tube Furnace

- Variable Speed Control
- Temperature Control System



## eXPRESS-LINE Solid Tube Furnace

- Available in 1" To 3" ID up to 24" Long
- Includes Touchscreen SmartControl